

# ABUS System Consumable Supplies

GE Healthcare offers an array of consumable supplies for the ABUS (Automated Breast Ultrasound) systems including scan membranes, POLYSONIC® ultrasound lotions and Protex™ disinfectants.

## ABUS Single-Use Stabilization Membrane for the Invenia™ ABUS and somo•v™ ABUS Platinum

The ABUS disposable scan membrane provides comfort to the patient during the 3D image acquisition. The flexible, polyester, mesh fabric stabilizes the breast tissue and enables better image quality during the scan. The membrane is single-use-only which eliminates the risk of cross-contamination and the need for sterilization.

Curved membrane dimensions: 229mm (L) x 213mm (W) x 19mm (D)

Catalog #	Description	Quantity per Box
E8340AB	ABUS Single-Use Stabilization Membrane for Invenia™ ABUS and somo•v™ ABUS Platinum models (non-Sterile)	50/box



## ABUS Single-Use Stabilization Membrane for the somo•v™ ABUS Insight

The ABUS disposable scan membrane provides comfort to the patient during the 3D image acquisition. The flexible, polyester, mesh fabric stabilizes the breast tissue and enables improved image quality during the scan. The membrane is single-use-only which eliminates the risk of cross-contamination and the need for sterilization.

Flat membrane dimensions: 229mm (L) x 213mm (W) x 19mm (D)

Catalog #	Description	Quantity per Box
E8340AC	ABUS Single-Use Stabilization Membrane for somo•v™ ABUS Insight models (non-Sterile)	50/box



## PolySonic Lotion

POLYSONIC® Ultrasound Lotion is an acoustic coupling lotion that is hypo allergenic and has a light pleasing scent.

Catalog #	Description
E8365EC	PolySonic Ultrasound Lotion (Non-Sterile) – 8.5 oz dispenser bottles
E8365EF	PolySonic Ultrasound Lotion – (Non-Sterile) 1 US gallon with refillable dispenser



## Protex Disinfectant Wipes

Protex™ one-step cleaner disinfectant kills more than 30 pathogens including HIV, H1N1 and MRSA, while remaining safe for even the most sensitive surfaces. The pre-moistened wipes are easy and convenient to use.

Catalog #	Description
E8365EG	Protex Ultra Disinfectant Wipes – 100 ct canister pack – 7 in x 10 in
E8365EB	Protex Ultra Disinfectant Wipes – 60 ct Soft-pack w/stay-put adhesive backing – 6.5in x 6.0in



## For additional information or to order contact:

GE Healthcare Accessories & Supplies

Phone: 866-281-7545

Fax: 877-279-6990

Email: [gehcaccessorysales@ge.com](mailto:gehcaccessorysales@ge.com)

[www.gehealthcare.com/accessories](http://www.gehealthcare.com/accessories)

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue any products described at any time without notice or obligation. Please contact your GE representative for the most current information.

© 2015 General Electric Company-All rights reserved.

GE, GE Monogram and imagination at work are trademarks of General Electric Company.

POLYSONIC and Protex are trademarks of Parker Laboratories, Inc.

\*Trademarks of General Electric Company.

GE Healthcare, a division of General Electric Company.

GE Healthcare  
3000 North Grandview Blvd  
Waukesha, WI 53188  
U.S.A

JB18299US(1)  
ACCD-0232-03.15-EN-US



imagination at work